

# Medilux Healthcare Ltd.

Healthcare for the Internet Age

[www.mediluxhealth.net](http://www.mediluxhealth.net)  
[www.mediluxprofessional.net](http://www.mediluxprofessional.net)

Safe - Natural - Drug-free - Non-invasive

## Q.Light® Phototherapy from Switzerland

The full range from Q.Products AG (Biotechnology & Photomedicine)

### Flexible Systems (add any of the filter modules below)

**Professional Users** – commercial premises  
 Medical: hospitals, clinics, nursing homes or salons  
 Cosmetic/colour: salons, treatment centres

#### Q.Light® PRO



Light diameter max 40cm. Aperture control, timer, function display

#### Q.Light® Skin Care



Light diameter max 25cm. On/off control

#### Filter modules

Q.Light® PRO and Q.Light® Skin Care require filter modules. Both accept all the filter modules shown.

#### Stand PRO




#### Stand Regular




### Specialised Systems (complete)

**Professional Users** - mobile  
 Rentals to patients  
**Home Use** - private purchase





**Photodynamic Therapy (PDT) for skin cancer**  
 Two filter modules: PDT Diagnosis and PDT Therapy  
 For medical professionals only – enquire for details



**Wound Care filters (2)**  
**With infra-red** for uninfected wounds, injuries, post-operative, muscle strains, cuts, grazes. **Without infra-red** for infected wounds, burns, pressure sores, leg ulcers, diabetic gangrene


**Q.Light® Wound Care**






**Pain Care filter**  
 For strains and injuries, back pain, rheumatism, arthritis, post operative pain, painful cosmetic procedures


**Q.Light® Pain Care**






**Acne Care filter**  
 For acne papulopustulosa and acne conglobata, common acne, acne vulgaris, acne comedonica,


**Q.Light® Acne Care**






**Psoriasis Care filter**  
 Eases psoriasis, neurodermatitis and atopic dermatitis. Relieves itching and improving skin tone


**Q.Light® Psoriasis Care**






**S.A.D. Care filter**  
 For Seasonal Affective Disorder and depressions. White light 420-780nm, free from harsh flicker of neon tubes in many 'light boxes'.


**Q.Light® SAD Care**






**Colour filter sets**  
**Regular:** red, orange, yellow, green, blue, violet  
**Advanced:** pink, cyan, turquoise, magenta, indigo


**Q.Light® Home Care**



Multi-purpose Q.Light® for a wide range of medical and cosmetic treatments in the home. No filters required, simple on/off control. Ideal for rentals to patients or for family medical use from acne to rheumatism and for personal skin and beauty treatments.





**Skin Care filter**  
 For daily skin care, elasticity and enhanced moisturisation



**Anti-Ageing filter**  
 For reducing fine lines and wrinkles, stimulates tired skin

**Stand Home**





**Cellulite Care filter**  
 For treatment of cellulite and weak connective tissue

For details download product brochures from [www.mediluxprofessional.net](http://www.mediluxprofessional.net)

# Medilux Healthcare Ltd.

Healthcare for the Internet Age

www.mediluxhealth.net  
www.mediluxprofessional.net

Safe - Natural - Drug-free - Non-invasive

## Q.Light® Phototherapy - Accessories

from Q.Products AG (Biotechnology & Photomedicine), Switzerland

### The Stands

#### Stand 'Home'



##### Accepts models

Acne Care  
Pain Care  
Wound Care  
Psoriasis Care  
Home Care  
Skin Care

##### Suitable for

Private homes  
Travelling

#### Stand 'Regular'



##### Accepts models

All models

##### Suitable for

Clinics  
Beauty salons  
Private homes

#### Stand 'PRO-Plus' \*\*



##### Accepts models

All models

##### Suitable for

Hospitals  
Nursing homes  
Clinics  
Beauty salons

\*\* Stand 'PRO-Plus' now has additional lower shelf and wire storage basket

Unless specified otherwise, Q.Light® PRO will be supplied with the 'Regular' stand and all other Q.Lights® will be supplied with 'Home' stands. Please check to ensure that you order the correct stand for your purposes. You may buy another stand later if you wish but this will be a separate purchase.

### The Filter Module Sets

for Q.Light® PRO and Q.Light® Skin Care only



<b>Medical Professional</b>	Acne	Pain	Wound (exc infra-red)	Wound (inc infra-red)	S.A.D.	Psoriasis
<b>Nursing Professional</b>		Pain	Wound (exc infra-red)	Wound (inc infra-red)	S.A.D.	
<b>Sports Professional</b>		Pain	Wound (exc infra-red)	Wound (inc infra-red)		
<b>Photodynamic Therapy</b>		PDT Diagnosis	PDT Therapy	[Restricted to skin cancer specialists]		
<b>Cosmetic Home</b>			Skin Care	Anti-Ageing	Cellulite	
<b>Cosmetic Professional</b>	Acne	Pain	Skin Care	Anti-Ageing	Cellulite	
<b>Colour Regular</b>	Red	Orange	Yellow	Green	Blue	Violet
<b>Colour Advanced</b>	Pink	Magenta	Turquoise	Cyan	Indigo	
[Not recommended for untrained users]						

Q.Lights® are EU certified medical products.

For more details visit: [www.mediluxhealth.net](http://www.mediluxhealth.net) or [www.mediluxprofessional.net](http://www.mediluxprofessional.net) Enquiries: [office@mediluxhealth.net](mailto:office@mediluxhealth.net)

NOTE: Medilux Healthcare Ltd. offers products on the basis of the Manufacturers' representations as to quality and efficacy. Any further guidance as to conditions which may benefit from use of the goods is based purely on our own opinions and beliefs regarding general application. We do not guarantee that goods will be suitable or effective for all purchasers and we do not advise on individual cases. We recommend that customers seek qualified medical advice before commencing treatment.